
# Anipak<sup>®</sup>

Smart Series 220

## UV ANTIBACTERIAL TOOTHBRUSH STERILIZER

Sterilize Against  
**99.9%**  
Bacteria- Viruses

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## About of product

### EN

Toothbrushes can provide an appropriate substrate for the growth of microbes, fungi and bacteria. ANIPAK toothbrush sterilizer has been produced with the utilization of the UVC ray in order to hold and antiseptics of the toothbrushes. This toothbrush holder is portable, and as soon as the toothbrush is placed in it, the UVC lamp turns on and then after 5 minutes radiation on the toothbrush this lamp is turned off. During this period, the air enters the toothbrush from the holes on the body of the toothbrush and dries the toothbrush, and with the radiation of the UVC ray simultaneously, 99% of the existing microorganisms on the toothbrush are disinfected and eliminated. 2 AAA batteries are used to supply the energy needed for the lamp of this toothbrush holder.

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## What is UV light?

### EN

The sun reflects a lot of ultraviolet rays in three UVA, UVB, and UVC spectra. Ultraviolet ray is electromagnetic energy which has a short wavelength and much energy and is invisible to the human eye, And this ray is in the electromagnetic spectrum, between X and visible light rays. Because of the ultraviolet absorption feature in the ozone layer of the Atmosphere, 99% of the ultraviolet ray that reaches to the ground is a type of spectrum of UVA with lower risk. Ultraviolet ray of UV-A type (315 ~ 400nm): This spectrum of radiation ray is apparently less risky and has the most spectrum of the ultraviolet radiation ray.

UV-B Radiation ray (280 ~ 315nm): This range of ultraviolet radiation ray is hazardous and causes the sunburn of the people's skin.

UV-C radiation ray (200 ~ 280nm): This spectrum of radiation ray is very destructive and causes severe sunburn of the people's skin. Scientific research has proven that there is the most effective frequency for eliminating the microorganisms in the UV-C spectrum of the radiation ray. This spectrum of radiation ray causes the destruction of DNA and thus kills 99.9% of germs, bacteria, and fungi. This spectrum of radiation ray can destroy DNA, resulting in the loss of 99-90% of germs, bacteria, and fungi.

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


## The "INTERTEK LABTEST" UV-C light wand Antibacterial Activity Test

NAME OF TEST ORGANISM (STRAIN NUMBER)	Klebsiella pneumoniae (ATCC 4352)	Pseudomonas aeruginosa (ATCC 9027)	MRSA (ATCC 33591)	Salmonella typhi (ATCC 6539)	Staphylococcus aureus (ATCC 6538)	Escherichia coli (ATCC 8739)
INITIAL COUNT (CFU/mL)	2.14*10 <sup>5</sup>	5.5*10 <sup>4</sup>	1.07*10 <sup>5</sup>	6.4*10 <sup>4</sup>	2.0*10 <sup>5</sup>	1.75*10 <sup>5</sup>
UV LIGHT EXPOSURE TIME	5 min	5 min	1 min	45 sec	1 min	5 min
THE NUMBER OF TEST ORGANISM AFTER UV LIGHT EXPOSURE	<10	<10	80	60	<10	<10
PERCENT REDUCTION OF TEST ORGANISM (%)	<b>99.99%</b>	99.98%	99.93%	<b>99.91%</b>	<b>99.99%</b>	<b>99.99%</b>

## Sterilization

### EN

**Maximum Radiation:** The stainless steel mirror inside this unit radiates the UVC ray to areas of the toothbrush that these beams do not reach them, and maximizes the disinfection capabilities of this device.

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### EN


**Automatic Timer:** This device has an internal timer, and as soon as the case of this device is closed, the lamp turns on and UVC ray radiation begins, and the lamp turns off after 5 minutes automatically.

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## UV Technology

### EN

Using the UVC technology is the best and most effective method for antiseptis and sterilization of the toothbrush in the form of a toothbrush holder. An ultraviolet light bulb with strong effectiveness is used in this device that eliminates 99.9% of bacteria, and germs on the surface of the toothbrush.

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
# Anipak®

Smart Series 220

UV ANTIBACTERIAL  
TOOTHBRUSH STERILIZER

## Portable

EN

This device is compact and its weight is only 140 grams. It works with 2 AAA 1.5V batteries and can be easily placed in a suitcase or bag.

This portable device is excellent for use in home, travel, school, sports club, holidays, workplace, and other items like these.

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## Slim design

### EN

The design of this device is such that it occupies a little space, and in addition to disinfection capability, this toothbrush holder is a suitable case for keeping the most toothbrushes.

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### Design

#### EN

The body of this lightweight device has been made from the plastic with ABS quality in 4 colors and this device is resistant against the fingerprint and scratches.

The special holes for air conditioning in this toothbrush holder reduce the excess moisture and the lid of this device prevents from entering the dust and bacteria in this device.

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
**EN**

The MODJE Industrial Group was established in 1995 in an area of 2500 square meters in the field of designing and manufacturing of building sanitary ware, and through its years of operation, it has been Succeeded to produce more than 60 new products with excellent quality for domestic supply and export to neighboring countries, utilizing its expert manpower and modern facilities.

This production complex has research and development, design, die making, production, plating, assembling and packaging.

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No.22, 5th East alley, Bist Meyery Lale  
st., Khoramdasht Industry Complex,  
Abali Road, Tehran, IRAN  
Tel: +9821 76212666  
Mobile: +98 9129172565  
fax: +9821 76217484

[Anipak.ir](http://Anipak.ir)